

Release Notes

DISCLAIMER

This software has not undergone extensive formal testing at this time. It will be included in the next version of the ProteinPilot™ 4.0 software, at which time it will have been subjected to such testing. AB SCIEX assumes no liability or responsibility for any loss howsoever caused in connection with using this software.

ProteinPilot™ 4.0 Patch for Windows 7 SP1(Rev.2)

This patch is distributed to users of the ProteinPilot™ 4.0 and 4.1 beta software on Windows 7 operation system when the application stops functioning after Service Pack 1 is installed.

Where to Get Help

Support

For support, contact AB SCIEX.

- Email: support@absciex.com
- Web: www.absciex.com

Prerequisites

This patch requires the ProteinPilot™ 4.0 software or ProteinPilot 4.1 beta software to be installed on Windows 7 operation system with Service Pack 1

Installing the Software

Double-click **ProteinPilot™ 4.0 Patch for Windows 7 SP1(Rev.2).exe** and then follow the instructions in the installation wizard.

To verify installation:

Win 7 64-bit System: navigate to 'Start > Control Panel > Programs and Features' window.

Win 7 32-bit System: navigate to 'Start > Control Panel > Programs > Programs and Features' window.

The list of currently installed programs should include **ProteinPilot™ 4.0 Patch for Windows 7 SP1(Rev.2)**

Updated files

The patch updates the following ProteinPilot™ software file in the **rscoree.dll** version 4.1.0.0 in the default location:

Win 7 64-bit System: "C:\Program Files (x86)\AB SCIEX\ProteinPilot\"

Win 7 32-bit System: "C:\Program Files\AB SCIEX\ProteinPilot\"

Uninstalling the Patch

1. Open the Control Panel and go to Programs and Features window (please see above)
2. Select the **ProteinPilot™ 4.0 Patch for Windows 7 SP1(Rev.2)** from the list and then click **Remove**.

Note: Uninstall will remove the installation record, but rscoree.dll is not removed from the installation directory.

Known Limitations

None

This document is provided to customers who have purchased AB SCIEX equipment to use in the operation of such AB SCIEX equipment. This document is copyright protected and any reproduction of this document or any part of this document is strictly prohibited, except as AB SCIEX may authorize in writing.

Software that may be described in this document is furnished under a license agreement. It is against the law to copy, modify, or distribute the software on any medium, except as specifically allowed in the license agreement. Furthermore, the license agreement may prohibit the software from being disassembled, reverse engineered, or decompiled for any purpose.

Portions of this document may make reference to other manufacturers and/or their products, which may contain parts whose names are registered as trademarks and/or function as trademarks of their respective owners. Any such usage is intended only to designate those manufacturers' products as supplied by AB SCIEX for incorporation into its equipment and does not imply any right and/or license to use or permit others to use such manufacturers' and/or their product names as trademarks.

AB SCIEX makes no warranties or representations as to the fitness of this equipment for any particular purpose and assumes no responsibility or contingent liability, including indirect or consequential damages, for any use to which the purchaser may put the equipment described herein, or for any adverse circumstances arising therefrom.

For Research Use Only. Not for use in diagnostic procedures.

The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners.

AB SCIEX™ is being used under license.

AB SCIEX

71 Four Valley Dr., Concord, Ontario, Canada. L4K 4V8.

AB SCIEX is ISO 9001 registered.

© 2011 AB SCIEX.

Printed in Canada

